

A

CONTROL III®
Disinfectant Germicide

Bulletins:

For
Home Care

And

Efficacy of
Control III® Products
For the Control of the
Human Immunodeficiency Virus Type 1
(AIDS Virus)

- Odorless
- Time Tested Since 1974
- EPA Registered Disinfectant

WHAT IS CONTROL III DISINFECTANT GERMICIDE?

This product is a time-proven germicide. It is an excellent disinfectant for inanimate hard surfaces and noncritical equipment in the home. This product, when used in an overall cleaning program, will help prevent cross contamination from *Pseudomonas aeruginosa*; *Staphylococcus aureus*; *Staphylococcus aureus* – MRSA; *Salmonella choleraesuis*, *Escherichia coli*, Herpes Simplex virus type 1, Herpes Simplex virus type e, Influenza virus type A, Vaccinia virus, Adenovirus type 5 and HIV-1. 1

HOW TO PREPARE CONTROL III® DISINFECTANT GERMICIDE

Prepare a disinfecting solution of CONTROL III DISINFECTANT GERMICIDE by adding ½ ounce (1 Tablespoon) of this product to 2 quarts of water in a plastic container.** This solution may be reused for up to 14 days. Solution may be checked daily with test strips to ensure adequate disinfection level of this solution. If the solution becomes visibly dirty or falls below 1,000 ppm (parts per million), dump the used solution and make a fresh solution. Keep disinfecting solution covered when not in use.

Bacterial Stability of Use Solution: Tests confirm that the use solution remains effective against the organisms listed on the label for up to 14 days when stored in a sealed container such as a spray bottle. If the product becomes visibly dirty or contaminated, the use-dilution must be discarded and a fresh product prepared. Always use clean, properly labeled containers when diluting this product. Bacterial stability of use solutions does not apply to open containers such as buckets or pails without lids. If using a bucket or pail without a lid, prepare a fresh solution daily or more often if the solution becomes visibly dirty or diluted.

**Re-use container must be labeled with the product name (Re-use solution of Control III Disinfectant Germicide), signal word, preparation date (mm/dd/yy) and expiration date (mm/dd/yy).

DISINFECTION OF SURFACES AND NONCRITICAL EQUIPMENT IN THE HOME WITH CONTROL III DISINFECTANT GERMICIDE

1. Disassemble equipment according to manufacturer's recommendations.
2. Wash all surfaces thoroughly in mild detergent (such as liquid dishwashing detergent) and warm tap water. DO NOT USE SOAP. Rinse parts with tap water.
3. Immerse equipment parts in solution of CONTROL III DISINFECTANT GERMICIDE, ensuring that there are no trapped air bubbles and all surfaces are wet. Soak for 10 minutes.
4. Rinse all equipment parts with hot tap water (allow water to run for 3 minutes before rinsing parts) and allow to air dry. After parts are dry, cover them with a clean towel until they are ready for use.

Control III Disinfectant Germicide kills the HIV-1 (AIDS Virus) on hard, non-porous surfaces in the presence of organic soil (5% blood serum) when diluted at a minimum level of 1,000 ppm active quaternary for a contact time of 10 minutes at room temperature.

¹ CONTROL III DISINFECTANT GERMICIDE has not been tested for effectiveness against *Mycobacterium tuberculosis* and must not be relied upon when a tuberculocidal product is desired.

Control III Disinfectant Germicide diluted at 1 ounce per gallon of water equals 1562 ppm active quaternary.

CONTROL OF THE AIDS VIRUS

KILLS HIV ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS

Personal Protection: Clean-up should always be done wearing protective latex gloves, masks, and eye protection.

Cleaning Procedure: Blood and other body fluids containing HIV must be thoroughly cleaned from surfaces and objects before application of this product.

Contact Time: Allow surface to remain wet for 10 minutes.

Disposal of Infectious Materials: Bloody body fluids, cleaning materials and clothing should be autoclaved and disposed of according to local regulations for infectious waste disposal.

Maril Products, Inc.

320 West 6th Street
Tustin, CA 92780

800-546-7711
714-544-7711

control3@earthlink.net

EPA No. 55364-3