




3-in-1 Clean, Disinfect, Sanitize

Handyclean STERIDOL® wipes are professional grade wipes with a powerful cleaning formula that includes unique active ingredients. The fungicidal, bactericidal, virucidal formula is EPA-approved, does not contain harsh solvents, and is proven effective against killing 99.9% of microorganisms tested. Each pre-moistened wipe delivers the correct amount of solution each time, quickly treating a variety of surfaces including, but not limited to: hospitals and healthcare settings, hotels and spas, airports and transportation hubs, gyms and fitness clubs, nursing homes and day cares, schools and cafeterias, grocery stores and supermarkets.

STERIDOL® WIPES CAN BE USED AGAINST 2019 NOVEL CORONAVIRUS

STERIDOL® wipes have demonstrated effectiveness against viruses similar to 2019 novel coronavirus (SARS-CoV-2) on hard non-porous surfaces. Therefore, this product can be used against 2019 novel coronavirus when used in accordance with the directions for use against Hepatitis A virus on hard, non-porous surfaces.

FEATURES

- 
Effective
 STERIDOL® wipes are effective against MRSA, Salmonella and flu virus.
- 
Safe on Surfaces
 Alcohol, bleach and solvent free wipes are gentle on sensitive equipment.
- 
Powerful & Fast-Acting
 Cleans in 5 seconds or less, Sanitizes in only 15 seconds, Disinfects in 4-10 minutes.

STERIDOL® WIPES

Product #	Description	# / Case
32031969	Handyclean Steridol Wipes Canister 75ct	6



STERIDOL[®] Efficacy Sheet

Bacteria	Clorox	Lysol	Sani-Cloth Plus	Steridol
Pseudomonas aeruginosa (Pseudomonas)	4 min	10 min	3 min	4 min
Salmonella enterica (Salmonella)	4 min	10 min	3 min	4 min
Staphylococcus aureus (Staph)	4 min	10 min	3 min	4 min
Acinetobacter baumannii	N/A	N/A	3 min	4 min
Burkholderia cepacia	N/A	N/A	N/A	
Campylobacter jejuni (Campylobacter)	4 min	N/A	3 min	4 min
Corynebacterium ammoniagenes	N/A	N/A	N/A	4 min
Enterococcus faecalis (Enterococcus)	N/A	N/A	N/A	4 min
Enterococcus faecalis - Vancomycin Resistant (VRE)	N/A	N/A	N/A	4 min
Escherichia coli (E. coli)	N/A	N/A	N/A	4 min
Escherichia coli O 157:H7	4 min	10 min	3 min	4 min
ESBL Escherichia coli - (Extended spectrum beta-lactamase producing E. coli)	N/A	N/A	N/A	4 min
Klebsiella pneumoniae (Klebsiella)	N/A	N/A	N/A	4 min
Legionella pneumophila	4 min	10 min	3 min	4 min
Listeria monocytogenes	N/A	N/A	N/A	4 min
Salmonella schottmuelleri (Salmonella)	4 min	N/A	3 min	4 min
Salmonella typhi (Salmonella)	N/A	N/A	N/A	4 min
Serratia marcescens (Serratia)	N/A	N/A	N/A	4 min
Shigella dysenteriae (Shigella)	N/A	N/A	N/A	4 min
Staphylococcus aureus - Community Associated Methicillin-Resistant (CA-MRSA) (NRS384) (USA30 0)	N/A	N/A	N/A	4 min
Staphylococcus aureus- Community Associated Methicillin-Resistant (CA-MRSA) (NRSI 23) (USA40 0)	N/A	N/A	N/A	4 min
Staphylococcus aureus- Methicillin-Resistant (MRSA)	N/A	N/A	3 min	4 min
Staphylococcus aureus- Multi-Drug Resistant (resistant to tetracycline Tc, penicillin Pc, streptomycin Sm) and erythromycin (Em) and susceptible to chloramphenicol (Cm) in vitro)	N/A	N/A	N/A	4 min
Staphylococcus aureus- Vancomycin Intermediate Resistant- (VISA); Streptococcus pyogenes (Strep)(a cause of scarlet fever)	N/A	N/A	N/A	4 min
Vibrio cholera	N/A	N/A	N/A	4 min
Viruses				
*Hepatitis B Virus (HBV)	N/A	N/A	3 min	4 min
*Hepatitis C Virus (HCV)	N/A	N/A	3 min	4 min
*Herpes Simplex Virus Type 1	4 min	10 min	3 min	4 min
*Herpes Simplex Virus Type 2	N/A	N/A	N/A	4 min
*HIV-1 (AIDS Virus)	4 min	10 min	3 min	1 min
*Human Coronavirus	4 min	10 min	3 min	4 min
*Inuenza Type A Virus / Brazil (Inuenza) ((the) Flu (Virus))	4 min	10 min	3 min	4 min
*Inuenza A Virus H 1 N 1	4 min	10 min	3 min	4 min
*2 013 Inuenza A Virus (H7N9)	N/A	N/A	N/A	4 min
*Norwalk Virus (Feline Calicivirus) (Norovirus)	N/A	N/A	N/A	10 min
*Respiratory Syncytial Virus (RSV)	4 min	10 min	3 min	4 min
*Rotavirus	4 min	N/A	N/A	10 min
*SARS Associated Coronavirus (SARS) (cause of Severe Acute Respiratory Syndrome)	N/A	N/A	N/A	10 min
*Vaccinia Virus (Pox Virus)	N/A	N/A	N/A	4 min
Animal Viruses				
Avian Inuenza Virus H3N2	N/A	N/A	N/A	4 min
Avian Inuenza Virus H5N1	4 min	10 min	3min	4 min
Canine Distemper Virus	N/A	N/A	N/A	4 min
Feline Calicivirus	N/A	N/A	N/A	10 min
Newcastle's Disease Virus	N/A	N/A	N/A	4 min
Pseudorabies Virus	4 min	N/A	N/A	4 min
Candida albicans	N/A	N/A	3 min	4 min
Mold & Mildew				
	N/A	N/A	N/A	4 min

The complete packaging label is available upon request. Information is based on what was readily available on respective web sites. Any and all trademarks are trademarks of respective owners. Information subject to change without notice.