

Data Sheet



Manufacturer:	eemagine Medical Imaging Solutions GmbH Gubener Str. 47 10243 Berlin Germany
Product Name:	EEG cap, Tin, 21 channels plus GND + REF.
Product Code:	CS-301
Long Description:	EEG cap, 21 electrodes plus GND + REF, reference position at CPz, layout 10/20, Tin electrode, 1x DB25 Male connector, unshielded cable 1.5m
Available sizes:	N1, N2, N3, N4, N5, B, I, C, S, M, L
Connector type:	DB 25
Connector gender:	Male
Cable length:	1.5 m

Quality Control Form

Serial nr:

Product testing passed:

Date:

Signature:

Pinning Scheme

Channel Number	Electrode	Pin numbers	
1	Fp1	1	
2	Fp2		14
3	F3	2	
4	F4		15
5	C3	3	
6	C4		16
7	P3	4	
8	P4		17
9	O1	5	
10	O2		18
11	F7	6	
12	F8		19
13	T3	7	
14	T4		20
15	T5	8	
16	T6		21
-	GND	9	
17	Cz		22
18	Fz	10	
19	Pz		23
20	M1	11	
21	M2		24
-	-	12	
-	-		25
22	REF	13	

Note: GND is located at AFz position. Reference is located at CPz.

Cap Layout

