



Installing your new animation pack for **Bioexplorer and Brainmaster**

1. Download the zip file “new_animation_pack-Bioexplorer-Brainmaster.zip”
2. Double click the zip file to open your default zip file extraction program



**Your extraction program may look different*

3. Click the “Extract” button at the top of the program window
4. Choose a folder to place the extracted “AVI” folder in



*Click here to choose a folder then click ok

5. To use the new animations in Bioexplorer or Brainmaster, navigate to the location you chose for the “AVI” folder when choosing a video for playback